

CON-LAB Conducted Immunity Lab Software

Product Overview

Ease of Use TDK Conducted Immunity Lab makes it simple to perform conducted immunity and bulk current injection (BCI) tests. TDK Conducted Immunity Lab can be configured to run conducted immunity tests according to test standards applicable to your products. The familiar graphical interface makes it easy to use without sacrificing flexibility or performance.

Performance TDK Conducted Immunity Lab enables you to test your products for immunity to conducted disturbances on their power, signal, and control lines, and to determine your products' immunity to the complex waveforms of induced bulk cable current.


Flexibility Enables you to create simple test files for others to run, or you can interact with the test process every step of the way.

Test Standards

- Automotive (BMW, Chrysler, ECE, European Auto Directive, Ford, GMW, Hyundai, ISO 11451, ISO 11452, JASSO, Nissan, Porsche, PSA, Renault, SAE)
- Aviation/Military (MIL-STD-461, RTCA DO-160)
- Commercial (IEC 61000-4-6, ETSI)
- Telecommunications (GSM, CDMA, WCDMA, UMTS, DECT, Cordless Phones, Bluetooth, WiFi)

Minimum System Requirements

- Windows XP, Vista, or 7
- USB port
- National Instruments GPIB controller



TDK Conducted Immunity Lab Software combines a powerful set of testing tools with an easy-to-use interface to reduce the complexity associated with immunity testing.

Data Table Features

- Configurable data display format
- Flexible configuration allowing any data parameters to be displayed in the table
- Export data to a text file
- Data tables may be copied and pasted into documents in other applications (e.g. Word, etc.)
- Add/delete/move columns
- Color coding of data values to indicate markers
- Sort data on any data parameter
- "Find" feature to locate data

CON-LAB Conducted Immunity Lab Software

Testing Features

- Real-time data display during tests
- Configurable test level and dwell time
- Supports a wide variety of signal generators, amplifiers, power meters, receivers, spectrum analyzers, base station simulators, switch modules, and other types of equipment
- E-mail notification of test completion with the test file optionally attached
- Configurable delays
- Analysis of measurement uncertainty per UKAS LAB34
- Tests may be paused, resumed, or stopped at any time
- Easy-to-edit factor and test level files
- Interactive control of the test to investigate the immunity threshold of the EUT
- Interacts with TDK EUT Monitoring software to monitor EUT parameters

Graph Features

- Display multiple graphs
- Flexible configuration allowing any data parameters to be plotted on any axis (x-axis, left y-axis, right y-axis)
- Linear or logarithmic axes
- Automatic or user-defined scaling of axes
- Configurable measurement units, trace thickness, color, and labels
- Graphs may be copied and pasted into documents in other applications (e.g. Word, etc.)
- Zoom, multiple zooms, nested zooms
- Add/delete color-coded data markers automatically or manually

Frequency Range Setup Features

- Frequency range may be split into sub-ranges to accommodate different equipment or settings in different frequency ranges
- Independent frequency stepping (linear, log, or file-based) in each sub-range
- Pause at the beginning or end of a sub-range and display a user-defined message
- Skip one or more sub-ranges

Ordering Information

Product: TDK Conducted Immunity Lab
Model Number: CON-LAB

To place an order or to learn more about the TDK products that best meet your needs, contact your TDK sales representative:

TDK RF Solutions Inc.

1101 Cypress Creek Rd.
Cedar Park, Texas 78613 USA
Phone: 1-512-258-9478
Fax: 1-512-258-0740
E-mail: info@tdkrf.com
World Wide Web: www.tdkrfsolutions.com

TOTAL RF EXPERTISE™



www.tdkrfsolutions.com

To learn more about TDK's wide range of innovative test products, solutions and services visit www.tdkrfsolutions.com

Copyright © TDK RF Solutions Inc. All rights reserved. Specifications subject to change without notice.

Rev. 2011-08-03